




Phospho-PAK1 (S204) Antibody

Product Code	CSB-PA020205
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13153
Immunogen	Synthesized peptide derived from Human PAK α around the phosphorylation site of S204.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PAK1; Serine/threonine-protein kinase PAK 1; Alpha-PAK; p21-activated kinase 1; PAK-1; p65-PAK
Immunogen Species	Homo sapiens (Human)
Target Names	PAK1

Image



Western Blot analysis of 3T3 cells using Phospho-PAK α (S204) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.