




Phospho-PARK2 (S131) Antibody

Product Code	CSB-PA020235
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60260
Immunogen	Synthesized peptide derived from Human Parkin around the phosphorylation site of S131.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PARK2; PRKN; E3 ubiquitin-protein ligase parkin; Parkinson juvenile disease protein 2; Parkinson disease protein 2
Immunogen Species	Homo sapiens (Human)
Target Names	PRKN

Image



Western blot analysis of 293T AD293 HeLa using p-Parkin (S131) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.