



TBX22 Antibody

Product Code	CSB-PA020251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y458
Immunogen	Synthesized peptide derived from the N-terminal region of Human TBX22.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TBX22; TBOX22; T-box transcription factor TBX22; T-box protein 22
Immunogen Species	Homo sapiens (Human)
Target Names	TBX22
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.