



ZDHHC2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA020268 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UIJ5 |
| Immunogen | Synthesized peptide derived from the Internal region of Human ZDHHC2. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat, Monkey |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:40000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | ZDHHC2; REAM; REC; ZNF372; Palmitoyltransferase ZDHHC2; Reduced expression associated with metastasis protein; Ream; Reduced expression in cancer protein; Rec; Zinc finger DHHC domain-containing protein 2; DHHC-2; Zinc finger protein 372 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ZDHHC2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |