



TCEB3B/TCEB3C Antibody

Product Code	CSB-PA020269
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYF1/Q8NG57
Immunogen	Synthesized peptide derived from the Internal region of Human Elongin A2/3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TCEB3B; TCEB3L; RNA polymerase II transcription factor SIII subunit A2; Elongin-A2; EloA2; Transcription elongation factor B polypeptide 3B; TCEB3C; TCEB3L2; RNA polymerase II transcription factor SIII subunit A3; Elongin-A3; EloA3; Transcr
Immunogen Species	Homo sapiens (Human)
Target Names	ELOA2/ELOA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.