



# PDGFRA Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA020307   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | P16234   |
| <b>Immunogen</b>           | Synthesized peptide derived from Human PDGFR- $\alpha$ around the non-phosphorylation site of Y754.  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human, Mouse, Rat  |
| <b>Tested Applications</b> | IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.  |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.  |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | PDGFRA; PDGFR2; RHEPDGFRA; Platelet-derived growth factor receptor alpha; PDGF-R-alpha; PDGFR-alpha; Alpha platelet-derived growth factor receptor; Alpha-type platelet-derived growth factor receptor; CD140 antigen-like family member A; CD14 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | PDGFRA   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.  |