



Phospho-HDAC7 (S155) Antibody

Product Code	CSB-PA030002
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUI4
Immunogen	Synthesized peptide derived from Human HDAC7 around the phosphorylation site of S155.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC7; HDAC7A; Histone deacetylase 7; HD7; Histone deacetylase 7A; HD7a
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.