



TMBIM4 Antibody

Product Code	CSB-PA030016
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC24
Immunogen	Synthesized peptide derived from the Internal region of Human GAAP.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TMBIM4; GAAP; LFG4; CGI-119; Protein lifeguard 4; Golgi anti-apoptotic protein; Protein S1R; Transmembrane BAX inhibitor motif-containing protein 4; Z-protein
Immunogen Species	Homo sapiens (Human)
Target Names	TMBIM4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.