



EMCN Antibody

Product Code	CSB-PA030025
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULC0
Immunogen	Synthesized peptide derived from the C-terminal region of Human Endomucin.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EMCN; EMCN2; MUC14; Endomucin; Endomucin-2; Gastric cancer antigen Ga34; Mucin-14; MUC-14
Immunogen Species	Homo sapiens (Human)
Target Names	EMCN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.