



ZAK Antibody

Product Code	CSB-PA030032
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYL2
Immunogen	Synthesized peptide derived from Human MLTK around the non-phosphorylation site of S165.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MLTK; ZAK; HCCS4; Mitogen-activated protein kinase kinase kinase MLT; Human cervical cancer suppressor gene 4 protein; HCCS-4; Leucine zipper and sterile alpha motif-containing kinase; MLK-like mitogen-activated protein triple kinase; Mixe
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.