



ATP5G3 Antibody

Product Code	CSB-PA030034
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48201
Immunogen	Synthesized peptide derived from the N-terminal region of Human ATP5G3.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ATP5G3; ATP synthase lipid-binding protein; mitochondrial; ATP synthase proteolipid P3; ATPase protein 9; ATPase subunit c
Immunogen Species	Homo sapiens (Human)
Target Names	ATP5G3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.