




Phospho-PGR (S294) Antibody

Product Code	CSB-PA030046
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06401
Immunogen	Synthesized peptide derived from Human PR around the phosphorylation site of S294.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
Immunogen Species	Homo sapiens (Human)
Target Names	PGR

Image



Western Blot analysis of 293 cells using Phospho-PR (S294) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.