



PHKA1/PHKA2 Antibody

Product Code	CSB-PA030051
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46020/P46019
Immunogen	Synthesized peptide derived from the N-terminal region of Human PHKA1/2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PHKA1; PHKA; Phosphorylase b kinase regulatory subunit alpha; skeletal muscle isoform; Phosphorylase kinase alpha M subunit; PHKA2; PHKLA; PYK; Phosphorylase b kinase regulatory subunit alpha, liver isoform; Phosphorylase kinase alpha L sub
Immunogen Species	Homo sapiens (Human)
Target Names	PHKA1/PHKA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.