



Phospho-PIM1 (Y309) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA030055 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P11309 |
| Immunogen | Synthesized peptide derived from Human Pim-1 around the phosphorylation site of Y309. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:40000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | PIM1; Serine/threonine-protein kinase pim-1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PIM1 |