



PIN1 Antibody

Product Code	CSB-PA030062
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13526
Immunogen	Synthesized peptide derived from Human Pin1 around the non-phosphorylation site of S16.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PIN1; Peptidyl-prolyl cis-trans isomerase NIMA-interacting 1; Peptidyl-prolyl cis-trans isomerase Pin1; PPlase Pin1; Rotamase Pin1
Immunogen Species	Homo sapiens (Human)
Target Names	PIN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.