



Phospho-PKD2 (S812) Antibody

Product Code	CSB-PA030064
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13563
Immunogen	Synthesized peptide derived from Human PKD2 around the phosphorylation site of S812.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PKD2; Polycystin-2; Autosomal dominant polycystic kidney disease type II protein; Polycystic kidney disease 2 protein; Polycystwin; R48321
Immunogen Species	Homo sapiens (Human)
Target Names	PKD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.