



Phospho-PLD1 (S561) Antibody

Product Code	CSB-PA030080
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13393
Immunogen	Synthesized peptide derived from Human PC-PLD1 around the phosphorylation site of S561.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLD1; Phospholipase D1; PLD 1; hPLD1; Choline phosphatase 1; Phosphatidylcholine-hydrolyzing phospholipase D1
Immunogen Species	Homo sapiens (Human)
Target Names	PLD1