




Phospho-YWHAQ (S232) Antibody

Product Code	CSB-PA030089
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27348
Immunogen	Synthesized peptide derived from Human 14-3-3 θ/τ around the phosphorylation site of S232.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	YWHAQ; 14-3-3 protein theta; 14-3-3 protein T-cell; 14-3-3 protein tau; Protein HS1
Immunogen Species	Homo sapiens (Human)
Target Names	YWHAQ

Image



Western blot analysis of 293T MOUSE-BRAIN using p-14-3-3 $\theta/$ (S232) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.