



P2RY13 Antibody

Product Code	CSB-PA030095
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPV8
Immunogen	Synthesized peptide derived from the Internal region of Human P2RY13.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	P2RY13; GPR86; GPR94; FKSG77; P2Y purinoceptor 13; P2Y13; G-protein coupled receptor 86; G-protein coupled receptor 94
Immunogen Species	Homo sapiens (Human)
Target Names	P2RY13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.