




FGFRL1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA030096 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8N441 |
| Immunogen | Synthesized peptide derived from the Internal region of Human FGFR-5. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:40000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | FGFRL1; FGFR5; FHFR; Fibroblast growth factor receptor-like 1; FGF receptor-like protein 1; FGF homologous factor receptor; FGFR-like protein; Fibroblast growth factor receptor 5; FGFR-5 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FGFRL1 |

Image



Western Blot analysis of 3T3 cells using FGFR-5 Polyclonal Antibody