



Phospho-CSNK1G1/CSNK1G2/CSNK1G3 (Y263) Antibody

Product Code	CSB-PA030100
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCP0/P78368/Q9Y6M4
Immunogen	Synthesized peptide derived from Human Casein Kinase I γ 1/2/3 around the phosphorylation site of Y263.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSNK1G1; Casein kinase I isoform gamma-1; CKI-gamma 1; CSNK1G2; CK1G2; Casein kinase I isoform gamma-2; CKI-gamma 2; CSNK1G3; Casein kinase I isoform gamma-3; CKI-gamma 3
Immunogen Species	Homo sapiens (Human)
Target Names	CSNK1G1/CSNK1G2/CSNK1G3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.