



PMS2CL/PMS2 Antibody

Product Code	CSB-PA030102
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68D20/P54278
Immunogen	Synthesized peptide derived from the C-terminal region of Human PMS2/PMS2CL.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PMS2CL; Protein PMS2CL; PMS2-C terminal-like protein; PMS2; PMSL2; Mismatch repair endonuclease PMS2; DNA mismatch repair protein PMS2; PMS1 protein homolog 2
Immunogen Species	Homo sapiens (Human)
Target Names	PMS2CL/PMS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.