



GPR85 Antibody

Product Code	CSB-PA030117
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60893
Immunogen	Synthesized peptide derived from the Internal region of Human GPR85.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR85; SREB2; Probable G-protein coupled receptor 85; Super conserved receptor expressed in brain 2
Immunogen Species	Homo sapiens (Human)
Target Names	GPR85