




ATR Antibody

Product Code	CSB-PA030128
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13535
Immunogen	Synthesized peptide derived from Human ATR around the non-phosphorylation site of S428.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ATR; FRP1; Serine/threonine-protein kinase ATR; Ataxia telangiectasia and Rad3-related protein; FRAP-related protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	ATR

Image



Western Blot analysis of COLO cells using ATR Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.