




APBB1IP Antibody

Product Code	CSB-PA030130
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5R6
Immunogen	Synthesized peptide derived from the Internal region of Human RIAM.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APBB1IP; PREL1; RARP1; RIAM; Amyloid beta A4 precursor protein-binding family B member 1-interacting protein; APBB1-interacting protein 1; Proline-rich EVH1 ligand 1; PREL-1; Proline-rich protein 73; Rap1-GTP-interacting adapter molecule; R
Immunogen Species	Homo sapiens (Human)
Target Names	APBB1IP

Image



Western Blot analysis of A549 cells using RIAM Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.