



POU4F1 Antibody

Product Code	CSB-PA030134
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01851
Immunogen	Synthesized peptide derived from the Internal region of Human Brn-3a.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	POU4F1; BRN3A; RDC1; POU domain; class 4, transcription factor 1; Brain-specific homeobox/POU domain protein 3A; Brain-3A; Brn-3A; Homeobox/POU domain protein RDC-1; Oct-T1
Immunogen Species	Homo sapiens (Human)
Target Names	POU4F1