



POU4F3 Antibody

Product Code	CSB-PA030135
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15319
Immunogen	Synthesized peptide derived from the C-terminal region of Human Brn-3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	POU4F3; BRN3C; POU domain; class 4, transcription factor 3; Brain-specific homeobox/POU domain protein 3C; Brain-3C; Brn-3C
Immunogen Species	Homo sapiens (Human)
Target Names	POU4F3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.