




PPARA Antibody

Product Code	CSB-PA030137
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07869
Immunogen	Synthesized peptide derived from Human PPAR- α around the non-phosphorylation site of S21.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPARA; NR1C1; PPAR; Peroxisome proliferator-activated receptor alpha; PPAR-alpha; Nuclear receptor subfamily 1 group C member 1
Immunogen Species	Homo sapiens (Human)
Target Names	PPARA

Image



Western Blot analysis of 3T3 cells using PPAR- α Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.