




BCAS3 Antibody

Product Code	CSB-PA030148
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6U6
Immunogen	Synthesized peptide derived from the Internal region of Human BCAS3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BCAS3; Breast carcinoma-amplified sequence 3; GAOB1
Immunogen Species	Homo sapiens (Human)
Target Names	BCAS3

Image



Western Blot analysis of 293 cells using BCAS3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.