




CDKAL1 Antibody

Product Code	CSB-PA030152
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VV42
Immunogen	Synthesized peptide derived from the N-terminal region of Human CDKAL1.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDKAL1; Threonylcarbamoyladenosine tRNA methyltransferase; CDK5 regulatory subunit-associated protein 1-like 1; tRNA-t(6)A37 methyltransferase
Immunogen Species	Homo sapiens (Human)
Target Names	CDKAL1

Image



Western Blot analysis of CoLo cells using CDKAL1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.