




SPDL1 Antibody

Product Code	CSB-PA030153
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EA4
Immunogen	Synthesized peptide derived from the C-terminal region of Human CCDC99.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCDC99; Protein Spindly; hSpindly; Arsenite-related gene 1 protein; Coiled-coil domain-containing protein 99; Rhabdomyosarcoma antigen MU-RMS-40.4A
Immunogen Species	Homo sapiens (Human)
Target Names	SPDL1

Image



Western Blot analysis of 3T3 cells using CCDC99 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.