



CSGALNACT2 Antibody

Product Code	CSB-PA030199
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6G5
Immunogen	Synthesized peptide derived from the Internal region of Human CsGalNAcT-2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSGALNACT2; CHGN2; GALNACT2; Chondroitin sulfate N-acetylgalactosaminyltransferase 2; Chondroitin beta-1; 4-N-acetylgalactosaminyltransferase 2; Beta4GalNAcT-2; GalNAcT-2
Immunogen Species	Homo sapiens (Human)
Target Names	CSGALNACT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.