



ARFGAP1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA030226 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8N6T3 |
| Immunogen | Synthesized peptide derived from the Internal region of Human ARF GAP1. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat, Monkey |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | ARFGAP1; ARF1GAP; ADP-ribosylation factor GTPase-activating protein 1; ARF GAP 1; ADP-ribosylation factor 1 GTPase-activating protein; ARF1 GAP; ARF1-directed GTPase-activating protein |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ARFGAP1 |