




PBK Antibody

Product Code	CSB-PA030268
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KB5
Immunogen	Synthesized peptide derived from Human PBK around the non-phosphorylation site of T9.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PBK; TOPK; Lymphokine-activated killer T-cell-originated protein kinase; Cancer/testis antigen 84; CT84; MAPKK-like protein kinase; Nori-3; PDZ-binding kinase; Spermatogenesis-related protein kinase; SPK; T-LAK cell-originated protein kinas
Immunogen Species	Homo sapiens (Human)
Target Names	PBK

Image



Western Blot analysis of 293T cells using PBK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.