





PRKG2 Antibody

Product Code	CSB-PA030289
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13237
Immunogen	Synthesized peptide derived from Human cGKII around the non-phosphorylation site of S126.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PRKG2; PRKGR2; cGMP-dependent protein kinase 2; cGK 2; cGK2; cGMP-dependent protein kinase II; cGKII
Immunogen Species	Homo sapiens (Human)
Target Names	PRKG2

Image



Western Blot analysis of A549 cells using cGKII Polyclonal Antibody



Western Blot analysis of 293T cells using cGKII Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.