




Phospho-MAP2K7 (S271) Antibody

Product Code	CSB-PA040000
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14733
Immunogen	Synthesized peptide derived from Human MEK-7 around the phosphorylation site of S271.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP2K7; JNKK2; MEK7; MKK7; PRKMK7; SKK4; Dual specificity mitogen-activated protein kinase kinase 7; MAP kinase kinase 7; MAPKK 7; JNK-activating kinase 2; MAPK/ERK kinase 7; MEK 7; Stress-activated protein kinase kinase 4; SAPK kinase 4; S
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K7

Image



Western Blot analysis of 3T3 cells using Phospho-MEK-7 (S271) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.