





SPHK2 Antibody

Product Code	CSB-PA040027
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRA0
Immunogen	Synthesized peptide derived from Human SphK2 around the non-phosphorylation site of T614.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SPHK2; Sphingosine kinase 2; SK 2; SPK 2
Immunogen Species	Homo sapiens (Human)
Target Names	SPHK2

Image



Western Blot analysis of 22RV1 cells using SphK2 Polyclonal Antibody



Western Blot analysis of AD293 using SphK2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.