




PAK7/PAK6 Antibody

Product Code	CSB-PA040042
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P286/Q9NQJ5
Immunogen	Synthesized peptide derived from Human PAK5/6 around the non-phosphorylation site of S602/S560.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PAK7; KIAA1264; PAK5; Serine/threonine-protein kinase PAK 7; p21-activated kinase 5; PAK-5; p21-activated kinase 7; PAK-7; PAK6; PAK5; Serine/threonine-protein kinase PAK 6; PAK-5; p21-activated kinase 6; PAK-6
Immunogen Species	Homo sapiens (Human)
Target Names	PAK5/PAK6

Image



Western Blot analysis of COLO cells using PAK5/6 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.