




ARID1B Antibody

Product Code	CSB-PA040080
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFD5
Immunogen	Synthesized peptide derived from the Internal region of Human BAF250b.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ARID1B; BAF250B; DAN15; KIAA1235; OSA2; AT-rich interactive domain-containing protein 1B; ARID domain-containing protein 1B; BRG1-associated factor 250b; BAF250B; BRG1-binding protein hELD/OSA1; Osa homolog 2; hOsa2; p250R
Immunogen Species	Homo sapiens (Human)
Target Names	ARID1B

Image



Western Blot analysis of 3T3 cells using BAF250b Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.