



Phospho-ARHGAP31 (T789) Antibody

Product Code	CSB-PA040084
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULL6
Immunogen	Synthesized peptide derived from Human CdGAP around the phosphorylation site of T789.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ARHGAP31; CDGAP; KIAA1204; Rho GTPase-activating protein 31; Cdc42 GTPase-activating protein
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGAP31
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.