




PEX5 Antibody

Product Code	CSB-PA040128
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50542
Immunogen	Synthesized peptide derived from the C-terminal region of Human Peroxin 5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PEX5; PXR1; Peroxisomal targeting signal 1 receptor; PTS1 receptor; PTS1R; PTS1-BP; Peroxin-5; Peroxisomal C-terminal targeting signal import receptor; Peroxisome receptor 1
Immunogen Species	Homo sapiens (Human)
Target Names	PEX5

Image



Western Blot analysis of A549 cells using Peroxin 5 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.