




RAB2A Antibody

Product Code	CSB-PA040134
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61019
Immunogen	Synthesized peptide derived from the Internal region of Human Rab 2A.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAB2A; RAB2; Ras-related protein Rab-2A
Immunogen Species	Homo sapiens (Human)
Target Names	RAB2A

Image



Western Blot analysis of A549 cells using Rab 2A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.