





RAB6A Antibody

Product Code	CSB-PA040135
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20340
Immunogen	Synthesized peptide derived from the Internal region of Human Rab 6A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAB6A; RAB6; Ras-related protein Rab-6A; Rab-6
Immunogen Species	Homo sapiens (Human)
Target Names	RAB6A

Image



Western Blot analysis of COLO cells using Rab 6A Polyclonal Antibody



Western Blot analysis of COLO NIH-3T3 cells using Rab 6A Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.