




Phospho-RAD51 (T309) Antibody

Product Code	CSB-PA040142
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06609
Immunogen	Synthesized peptide derived from Human Rad51 around the phosphorylation site of T309.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAD51; RAD51A; RECA; DNA repair protein RAD51 homolog 1; HsRAD51; hRAD51; RAD51 homolog A
Immunogen Species	Homo sapiens (Human)
Target Names	RAD51

Image



Western blot analysis of HELA-UV 293T-UV using p-Rad51 (T309) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.