



Phospho-RB1 (S249) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA040171 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P06400 |
| Immunogen | Synthesized peptide derived from Human Rb around the phosphorylation site of S249. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | IF, ELISA; IF:1:200-1:1000, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RB1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |