



RBBP8 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA040181 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q99708 |
| Immunogen | Synthesized peptide derived from Human CtIP around the non-phosphorylation site of S327. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | RBBP8; CTIP; DNA endonuclease RBBP8; CtBP-interacting protein; CtIP; Retinoblastoma-binding protein 8; RBBP-8; Retinoblastoma-interacting protein and myosin-like; RIM; Sporulation in the absence of SPO11 protein 2 homolog; SAE2 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RBBP8 |