




HRH4 Antibody

Product Code	CSB-PA040185
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3N8
Immunogen	Synthesized peptide derived from the Internal region of Human Histamine H4 Receptor.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HRH4; GPCR105; Histamine H4 receptor; H4R; HH4R; AXOR35; G-protein coupled receptor 105; GPRv53; Pfi-013; SP9144
Immunogen Species	Homo sapiens (Human)
Target Names	HRH4

Image



Western Blot analysis of COLO cells using Histamine H4 Receptor Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.