



Phospho-BCL2L1 (T115) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA040207 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q07817 |
| Immunogen | Synthesized peptide derived from Human Bcl-x around the phosphorylation site of T115. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | BCL2L1; BCL2L; BCLX; Bcl-2-like protein 1; Bcl2-L-1; Apoptosis regulator Bcl-X |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | BCL2L1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |