



REV3L Antibody

Product Code	CSB-PA040209
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60673
Immunogen	Synthesized peptide derived from the N-terminal region of Human DNA pol ζ.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	REV3L; POLZ; REV3; DNA polymerase zeta catalytic subunit; Protein reversionless 3-like; REV3-like; hREV3
Immunogen Species	Homo sapiens (Human)
Target Names	REV3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.