



BAZ2A Antibody

Product Code	CSB-PA040227
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIF9
Immunogen	Synthesized peptide derived from the Internal region of Human Tip5.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BAZ2A; KIAA0314; TIP5; Bromodomain adjacent to zinc finger domain protein 2A; Transcription termination factor I-interacting protein 5; TTF-I-interacting protein 5; Tip5; hWALp3
Immunogen Species	Homo sapiens (Human)
Target Names	BAZ2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.